

Perfecting Performance



Occlutech Stiff Guide Wire




The obvious leader

The Occlutech Stiff Guide Wire has been developed using a stainless-steel core and a unique PTFE coating. This combination gives well-balanced support for day-to-day cases, as well as advanced ones.

A soft, atraumatic J-tip, in combination with the exceptional torque control, allows exceptional navigation in difficult and challenging geometry.

This stiff configuration coupled with the flexible tip allows a good push while making the wire very safe and easy to use.

- **Soft atraumatic J-tip** – allowing safe navigation
- **Stiff configuration** – supporting easy entry into the left atrium
- **Completely coated** – minimizing the risk of thrombus in the left atrium



Product specifications

Article No.	Ø [inch]	Length [cm]	Tip Shape [mm]
52GWR01	0.035	260	J - 3 mm

