




Occlutech ASD Occluder

Compatibility chart & ordering information




The Procedure Pack
INCLUDES
The Accessory Medical Devices
Occlutech Pistol Pusher



Ø LA DISC: LEFT ATRIAL DISC DIAMETER
Ø WAIST: MAX. WAIST DIAMETER
Ø RA DISC: RIGHT ATRIAL DISC DIAMETER

Procedure Pack Art. No.	AS Defect Size (D) [mm]	Ø Waist [mm]	Ø LA Disc [mm]	Ø RA Disc [mm]	ODS* Size [F]	Occlutech Pistol Pusher (OPP) Color
29ASD04P	D≤4	4	11	9	7	●
29ASD05P	4<D≤5	5	14	11	7	●
29ASD06P	5<D≤6	6	16.5	12.5	7	●
29ASD07P	6<D≤7.5	7.5	18	14	7	●
29ASD09P	7.5<D≤9	9	20.5	16.5	7	●
29ASD10P	9<D≤10.5	10.5	22	18	7	●
29ASD12P	10.5<D≤12	12	27	23	9	●
29ASD13P	12<D≤13.5	13.5	28.5	24.5	9	●
29ASD15P	12<D≤15	15	30	26	9	●
29ASD16P	15<D≤16.5	16.5	31.5	27.5	9	●
29ASD18P	15<D≤18	18	33	29	9	●
29ASD19P	16.5<D≤19.5	19.5	34.5	30.5	10	●
29ASD21P	18<D≤21	21	36	32	11	●
29ASD24P	21<D≤24	24	39	35	11	●
29ASD27P	24<D≤27	27	42	38	12	●
29ASD30P	27<D≤30	30	45	41	12	●
29ASD33P	30<D≤33	33	48	43	12	●
29ASD36P	33<D≤36	36	52	46	12	●
29ASD39P	36<D≤39	39	54	49	12	●
29ASD40P	39<D≤40	40	55	50	12	●

Availability subject to regulatory approval. Product data may be subject to change. Contact your Occlutech representative for information.
* The Occluders are compatible with Occlutech Delivery Set, e.g. ODS 45° (Art. nN. 51DSXXX). Recommended size. Sold separately.



Occlutech International AB
Landskronavägen 2, SE-252 32 Helsingborg, Sweden
order@occlutech.com, www.occlutech.com

