



Occlutech PDA Occluder

Compatibility chart & ordering information





The Procedure Pack
INCLUDES
The Accessory Medical Devices
Occlutech Occlusion Pusher



D1: DISC DIAMETER
D2: WAIST DIAMETER
D3: SHANK DIAMETER

L: DEVICE LENGTH
AO: AORTA
PA: PULMONARY ARTERY

Procedure Pack Art. No.	Length [mm]	Ø D1 [mm]	Ø D2 [mm]	Ø D3 [mm]	ODS* Size [F]	Occlutech Occlusion Pusher (OP)
Standard Shank 						
42PDA050	4.25	9	3.5	5	6	500P120
42PDA060	5.00	10	4	6	6	500P120
42PDA070	6.05	11	5	7	6	500P120
42PDA080	6.30	13	6	8	6	500P120
42PDA100	7.00	16	8	10	7	500P020
42PDA120	12.00	18	10	12	7	500P020
42PDA150	14.00	20	12	15	8	500P020
42PDA180	16.00	24	14	18	9	500P020
Long Shank 						
43PDA05LO	7.00	9	3.5	5	6	500P120
43PDA06LO	7.50	10	4	6	6	500P120
43PDA07LO	8.50	11	5	7	6	500P120
43PDA08LO	9.00	13	6	8	6	500P120
43PDA10LO	10.50	16	8	10	7	500P020

Availability subject to regulatory approval. Product data may be subject to change. Contact your Occlutech representative for information.
* Occlutech Delivery Set (ODS). Include: ODS Version I - 180° - 80 cm.



Occlutech International AB
Landskronavägen 2, SE-252 32 Helsingborg, Sweden
order@occlutech.com, www.occlutech.com

