




Occlutech PLD Square

Compatibility chart & ordering information



The Procedure Pack
INCLUDES
The Accessory Medical Devices
Occlutech Flex Pusher I



- D1: PROXIMAL DISC LENGTH
- D2: DISTAL DISC LENGTH
- D3: WAIST DIAMETER (ONLY FOR SQUARE W*)
- D4: DIAGONAL LENGTH OF PROXIMAL DISC
- D5: DIAGONAL LENGTH OF DISTAL DISC
- H: HEIGHT OF THE WAIST (3 MM FOR ALL SIZES)

Procedure Pack Art. No.	[mm]	[mm]	[mm]	[mm]	[mm]	Occlutech Flex Pusher I (FP) Color	Introducing System ¹			
							Cook-Cook ² Loader	Sheath	Cook-Terumo ² Loader	Sheath
Square W*	D1	D2	D3*	D4	D5					
60PLD04WF	7.5	9	4	11.5	13	●	RCF-6F	KCFW-7F	RCF-7F	RSC04-7F
60PLD05WF	8.5	10	5	12.5	14	●	RCF-6F	KCFW-7F	RCF-7F	RSC04-7F
60PLD06WF	9.5	11	6	14	16	●	RCF-6F	KCFW-7F	RCF-7F	RSC04-7F
60PLD07WF	10.5	12	7	16	17	●	RCF-7F	KCFW-8F	RCF-7F	RSC04-7F
Square T*	Def. S	D1	D2	D4	D5					
62PLD03TF	3	7	8	10	11.5	●	RCF-6F	KCFW-7F	RCF-7F	RSC04-7F
62PLD05TF	5	8.5	10	12.5	14	●	RCF-6F	KCFW-7F	RCF-7F	RSC04-7F
62PLD07TF	7	11	12	16	17	●	RCF-7F	KCFW-8F	RCF-7F	RSC04-7F

Def. S = Defect Size

Availability subject to regulatory approval. Product data may be subject to change. Contact your Occlutech representative for information.

* W and T represent the types of connections between the discs. W stands for waist and T stands for twist.

1. Ordering Reference Numbers:

Loader:	French Size	Brand	Order No.	Ref No.	Length [cm]
	6F	Cook Performer (RCF)	G07337	RCF - 6.0 - 38 - J	13 cm
	7F	Cook Performer (RCF)	G07338	RCF - 7.0 - 38 - J	13 cm
	9F	Cook Performer (RCF)	G07564	RCF - 9.0 - 38 - J	13 cm
	10F	Cook Performer (RCF)	G07417	RCF - 10.0 - 38 - J	13 cm
Delivery Sheath:	French Size	Brand	Order No.	Ref No.	Length [cm]
	7F	Cook Ansel Flexor (KCFW)	G49044	KCFW - 7.0 - 35 - 110 - RB - HFANL1 - HC	110 cm
	7F	Terumo Destination (RSC04)	RSC04		90 cm
	8F	Cook Ansel Flexor (KCFW)	G49046	KCFW - 8.0 - 35 - 110 - RB - HFANL1 - HC	110 cm
	8F	Terumo Destination (54-89001)	54-89001		90 cm
	9F	Cook Performer (RCFW)	G09190	RCFW - 9.0 - 38 - 75 - RB - MTS	75 cm
	9F	Cook Performer (RCFW)	G08694	RCFW - 9.0 - 38 - 85 - RB - MTS	85 cm
	10F	Cook Performer (RCFW)	G08709	RCFW - 10.0 - 38 - 75 - RB - MTS	75 cm
	10F	Cook Performer (RCFW)	G08903	RCFW - 10.0 - 38 - 85 - RB - MTS	85 cm

2. Please note that Occlutech PLDs should not be retracted more than 3 times from the long sheath once they are deployed at the implantation site.




Occlutech International AB
Landskronavägen 2, SE-252 32 Helsingborg, Sweden
order@occlutech.com, www.occlutech.com




Occlutech PLD Rectangular



Compatibility chart & ordering information



The Procedure Pack
INCLUDES
The Accessory Medical Devices
Occlutech Flex Pusher I



AXB: WAIST LENGTH TO WIDTH (ONLY FOR RECTANGULAR W*)
D1XD3: PROXIMAL DISC LENGTH TO WIDTH
D2XD4: DISTAL DISC LENGTH TO WIDTH
D5: DIAGONAL LENGTH OF PROXIMAL DISC
D6: DIAGONAL LENGTH OF DISTAL DISC
H: HEIGHT OF THE WAIST (3 MM FOR ALL SIZES)

Procedure Pack Art. No.	[mm]	[mm]	[mm]	[mm]	[mm]	Occlutech Flex Pusher I (FP) Color	Introducing System ¹				
							Cook-Cook ² Loader	Sheath	Cook-Terumo ² Loader	Sheath	
 Rectangular W*	AxB	D1xD3	D2xD4	D5	D6						
61PLD04WF	4x2	8x6	9x7	10	11.5	●	RCF-6F	KCFW-7F	RCF-7F	RSC04-7F	
61PLD06WF	6x3	10x7	11x8	12.5	14	●	RCF-6F	KCFW-7F	RCF-7F	RSC04-7F	
61PLD08WF	8x4	12x8	14x10	15	16.5	●	RCF-7F	KCFW-8F	RCF-7F	RSC04-7F	
61PLD10WF	10x4	14x8	16x10	17	19	●	RCF-7F	KCFW-8F	RCF-7F	54-89001-8F	
61PLD12WF	12x5	16x9	17x10	19	21	●	RCF-9F	RCFW-9F	N/A	N/A	
61PLD14WF	14x6	19x11	21x13	22	24	●	RCF-9F	RCFW-9F	N/A	N/A	
61PLD16WF	16x8	20x12	23x14	24.5	26.5	●	RCF-10F	RCFW-10F	N/A	N/A	
61PLD18WF	18x10	22x14	24x16	26.5	28.5	●	RCF-10F	RCFW-10F	N/A	N/A	
 Rectangular T*	Def. S	D1xD3	D2xD4	D5	D6						
63PLD05TF	5x3	9x6	10x8	11.5	13	●	RCF-6F	KCFW-7F	RCF-7F	RSC04-7F	
63PLD07TF	7x4	11x8	12x9	14	16	●	RCF-6F	KCFW-7F	RCF-7F	RSC04-7F	
63PLD10TF	10x4	14x8	16x10	17	19	●	RCF-7F	KCFW-8F	RCF-7F	54-89001-8F	
63PLD12TF	12x5	16x9	17x11	19	21	●	RCF-9F	RCFW-9F	N/A	N/A	

Def. S = Defect Size

Note: For the 9F and 10F size Paravalvular Leak Devices there is no compatible Terumo Long Sheath.

Availability subject to regulatory approval. Product data may be subject to change. Contact your Occlutech representative for information.

* W and T represent the types of connections between the discs. W stands for waist and T stands for twist.

1. Ordering Reference Numbers:

Loader:	French Size	Brand	Order No.	Ref No.	Length [cm]
	6F	Cook Performer (RCF)	G07337	RCF - 6.0 - 38 - J	13 cm
	7F	Cook Performer (RCF)	G07338	RCF - 7.0 - 38 - J	13 cm
	9F	Cook Performer (RCF)	G07564	RCF - 9.0 - 38 - J	13 cm
	10F	Cook Performer (RCF)	G07417	RCF - 10.0 - 38 - J	13 cm
Delivery Sheath:	French Size	Brand	Order No.	Ref No.	Length [cm]
	7F	Cook Ansel Flexor (KCFW)	G49044	KCFW - 7.0 - 35 - 110 - RB - HFANL1 - HC	110 cm
	7F	Terumo Destination (RSC04)	RSC04		90 cm
	8F	Cook Ansel Flexor (KCFW)	G49046	KCFW - 8.0 - 35 - 110 - RB - HFANL1 - HC	110 cm
	8F	Terumo Destination (54-89001)	54-89001		90 cm
	9F	Cook Performer (RCFW)	G09190	RCFW - 9.0 - 38 - 75 - RB - MTS	75 cm
	9F	Cook Performer (RCFW)	G08694	RCFW - 9.0 - 38 - 85 - RB - MTS	85 cm
	10F	Cook Performer (RCFW)	G08709	RCFW - 10.0 - 38 - 75 - RB - MTS	75 cm
	10F	Cook Performer (RCFW)	G08903	RCFW - 10.0 - 38 - 85 - RB - MTS	85 cm

2. Please note that Occlutech PLDs should not be retracted more than 3 times from the long sheath once they are deployed at the implantation site.



Occlutech International AB
Landskronavägen 2, SE-252 32 Helsingborg, Sweden
order@occlutech.com, www.occlutech.com

